

These tables have been specially designed in order to minimise the problems associated with Repetitive Strain Injury. The multi position memory of the height control enables both taller and shorter people to use the table with comfort, both sitting or standing. All that is required is a simple memory set procedure for each individual or application using the table. The memorised heights can then be used over and over again. Alternatively, there is a continuously variable setting giving all heights between the minimum and maximum heights available.



Splash-proof height controller



Heavy Duty 4 Leg Table for Heavy Items Such as Safety Cabinets

These adjustable height tables also incorporate true parallel drive that ensures that even if the table is obliquely loaded the table top will always remain level and true during use.

Mountings for fan / filter modules, power points, equipment or storage racks may be fitted to the worktop as required. These ready-assembled tables can be supplied manufactured in either grade 304 or 316 stainless steel. Actuators are finished in painted steel but can be clad in stainless steel at extra cost.

Four legged versions are available with added loading capacity and even extra rigidity - heavy duty and ideal for supporting safety cabinets or LAF units, etc.

Please call us with your requirements - standard and custom models.

Typical Specification:

All-welded 240-grit polished 304 stainless steel tabletop with matching box stainless support frame with a parallel telescopic height adjusters generally to the literature. The base is fitted with plastic feet or casters. Where fitted, 2 casters are lockable.

<<< >>>

Below the tabletop there is a controller mounted, which can programme four set positions or allow height adjustment over the full travel.

Height range:- static models: 675 to 1125 mm mobile models: 745 to 1195 mm

Maximum lift for 2-leg model:- 200KG

<<< >>>

Each table is supplied with earthing leads, a 4-m mains cable with a 13A style plug and a hook below the table for storage.



Mobile trolleys



Downflow Pods



Downflow Benches & Storage

As part of our policy of continuous product improvement, the specifications are subject to change without notice.

Adjustable Height Tables

Ringstones Industrial Estate, Whaley Bridge, High Peak, SK23 7PD, UK
Tel: 01663 735494 Fax: 01663 732744 Web: www.qedscientific.co.uk

