



[Model GP50WM](#)

## MEDIA WARMING CABINETS

These are designed for use in X-Ray & CT theatres for safely warming bottles of contrast media to body temperature and also for warming ultrasound gel, etc. It is considered best practice to warm fluids prior to use to ensure patient comfort.

### DESCRIPTION

These cabinets will give long & trouble-free life in normal use. They are fitted with a digital temperature controller as well as an over-temperature safety cutout for product protection.

A flat perforated shelf ensures free air movement around the bottles to guarantee stable control of the internal temperature at 37c.

### CONSTRUCTION

The outer cabinets are made from sheet steel finished in easy-clean grey stove powder paint; stainless steel exterior cabinets are an option. The stainless steel chambers are more hygienic and durable, and offer excellent protection from corrosion. The internal design ensures optimum control and warm-up. An internal fan is fitted to all of these digitally controlled models.

### CONTROLS

All models are fitted with a digital PID temperature controller with the temperature fixed at a maximum of 37C. On request, other fixed or variable temperatures are available. An over temperature safety thermostat with warning lamp is always fitted.



[Model GP40WM/2S with two shelves](#)

### CONTRAST MEDIA WARMING CABINET RANGE

Model Ref.		Bottles approx	No. Of Shelves	External (mm) Ht x W x D	Internal (mm) Ht x W x D	Weight Kg*
GP30WM	Digital control	10	1	570 x 415 x 422	220 x 315 x 315	20
GP40WM	Digital control	20	1	670 x 415 x 422	320 x 315 x 315	25
GP40WM/2S	Digital control	40	2	670 x 415 x 422	320 x 315 x 315	26
GP50WM	Digital control	30	1	615 x 590 x 422	240 x 490 x 315	30
GP75WM	Digital control	70	1	623 x 590 x 567	248 x 490 x 460	35

Custom build sizes available

Add 50mm to external width to allow for door opening.

\* Empty

#### Options

All cabinets can be factory fitted with wall fixing brackets and audible alarm.

Other accessories include stands, either with or without lockable wheels.

Larger warming cabinets are also available. Please call for details of 'QED fluid warming cabinets'.

#### Other equipment available includes:

Fluid / lotion / blanket warming cabinets

Heated patient warming blankets

In-line fluid and blood warming devices