

All models have adjustable temperature control up to approx. +50°C. The cabinets are made from sheet steel finished in easy clean grey stove paint, and the wire shelves are adjustable to accommodate glassware and other items of different sizes. Wig drying models also available.

The two smaller models can be bench or wall sited. Models DC/113 & DC/226 have sliding glass doors and can be wall-mounted if space is at a premium.



Models DC/1000 (two door) and DC/540 (single door) are floor-standing with hinged toughened glass doors. They can be fitted with an interior fan and dial thermostat control for faster warming.

**WIG DRYING** cabinets are available as versions of the 540 & 1000 litre models. These have wheels, special controllers and safety cut outs fitted.

### MODEL RANGE

Model Ref.	Capacity Litres	Shelves Supplied	External (mm)			Internal (mm)		
			Ht	x W	x D	Ht	x W	x D
DC/113	113	2	545	x 730	x 360	460	x 728	x 341
DC/226	226	3	652	x 928	x 449	567	x 926	x 430
DC/540	540	3	1650	x 583	x 610	1550	x 581	x 600
DC/1000	1000	4	1650	x 1110	x 585	1550	x 1108	x 575

#### Fan Circulation (F)

This gives a stable, even temperature within the chamber. This option is ideal for drying products that are temperature sensitive (for example biological / pharmaceutical preparations and plastics).

#### Fan Extract (EX)

This is ideal when very wet glassware or products with high water content are to be dried.

#### Circulation with Extract (F/EX)

Recommended when drying large volumes of very wet, temperature-sensitive product.

Model	Standard model type	Model	Standard model type
DC/113	Convection	DC/1000	Convection
DC/226	Convection	DC/1000/F	Circulation fan
DC/540	Convection	DC/1000/EX	Extract fan
DC/540/F	Circulation fan	DC/1000/F/EX	Circulation & extract fans
DC/540/EX	Extract fan	<i>Also pass-through, stainless steel, spark-proof &amp; nitrogen purged versions of DC540 &amp; DC1000</i>	
DC/540/F/EX	Circulation + extract fans		

### WARRANTY

We are so confident in the overall reliability of our drying cabinets that we give a full 12 months parts and labour on-site warranty (UK mainland).

**OTHER PRODUCTS** in the QED Scientific range include incubators, ovens and laboratory refrigerators.

As part of our policy of continuous product improvement, the specifications are subject to change without notice