

Mod.	270	370	440	530	625
Capacità (litri) Capacity (litres)	270	370	440	530	625
Dimensioni (L x P x H cm) Dimensions (W x D x H cm)	60 x 65 x 164	60 x 60 x 192	60 x 70 x 192	72 x 70 x 207	72 x 80 x 207
Finitura esterna External finishing	W / S	W / S	W / S	W / S	W / S
Finitura interna Internal finishing	W / S	W / S	W / S	W / S	W / S
Isolamento (mm) Insulation (mm)	50	50	50	60	60
Rullini Rollers	4				
Piedi Feet		4	4	4	4
Porta isolata Insulated door	1	1	1	1	1
Illuminazione interna Internal lighting	•	•	•	•	•
Ripiani Shelves	2	3	3	3	3
Termostato elettronico con allarme Thermostat with alarm	•	•	•	•	•
Pannello di controllo PRO Control panel PRO	• PRO	• PRO	• PRO	• PRO	• PRO
Refrigerazione ventilata Forced air refrigeration	•	•	•	•	•
Sbrinatorio automatico Automatic defrosting	•	•	•	•	•
Refrigerante Refrigerant	ISCEON 89	ISCEON 89	ISCEON 89	ISCEON 89	ISCEON 89
Voltaggio (V/ph/Hz) Voltage (V/ph/Hz)	230/1/50	230/1/50	230/1/50	230/1/50	230/1/50

W = Bianco / White / Blanc / Weiß S = Acciaio inox / Stainless steel / Acier inox / Edelstahl

UK Distributor:
QED Scientific Ltd.

Ringstones Industrial Estate, Whaley Bridge, High Peak, SK23 7PD
Tel: 01663 735494 Fax: 01663 732744 Web; www.qedscientific.co.uk



-40C Super Freezers

EVERmed
MEDICAL REFRIGERATION

EVERmed s.r.l

Via Aldo Moro n°25 46020 Pegognaga (MN) - Italy

Telef. +39-0376-550828 Telefax +39-0376-550831 <http://www.evermed.it> mail:evermed@evermed.it

PDF Models -30°-40°

EVERmed
MEDICAL REFRIGERATION

*La ditta fornitrice si riserva il diritto di apportare ai propri prodotti quelle modifiche tecniche e stilistiche che riterrà necessarie o utili senza pregiudicare le caratteristiche essenziali. *Should it be necessary or useful without undermining any model essential features, technical and stylistic characteristics are subject to be changed without manufacturer's previous notice. *La Maison constructrice se réserve le droit d'opérer sans préavis toutes modifications techniques et stylistiques qu'elle considèrera nécessaires ou utiles sans compromettre les caractéristiques essentielles des produits. *Der Hersteller behält sich das Recht vor, alle technischen und das Aussehen betreffenden Änderungen vorzunehmen, die er für nützlich oder nötig hält, ohne die wesentlichen Merkmale der Erzeugnisse zu ändern.



PDF 625 S



PDF 625 W



The **standard controller** has a digital read-out . Programming of the internal temperature is via this controller and has bright red digits and 4 switches to allow control and regulation of all the working functions. Included is a visual/ acoustic alarm to alert users of deviations from the set temperature and also a test function to detect defective sensors.



The **PRO control panel** is an innovative concept of refrigerator control. It has a large back-lit LCD display with black digits and constantly monitors the internal temperature with an accuracy of 0.5C as well as the full working status of the fridge (compressor running, defrost, etc). 4 soft-touch buttons give access to all the various setting menus so that you can have complete control. The alarm system signals 8 different conditions (high and low temperature, door ajar, power failure, battery discharged, dirty condenser and sensors failure) thereby giving maximum security to the stored product. The alarm has a battery back-up so that the user is alerted even in a power failure situation. Total product security is assured with key password protection (with automatic activation) to avoid tampering and the alarm can be remotely signalled by the NO/NC contacts or RS485 port (optional software required). Alarm conditions can be muted with automatic re-activation. An internal memory for events and alarms guarantees full traceability of the storage conditions.



PDF 270 S



PDF 370 S



PDF 440 S



PDF 530 S



PDF 625 S

The **-40C PDF super-freezer series** is available in 2 different finishing materials - steel epoxy coated white colour or stainless steel. Models have capacities from 270 to 625 litres. The two levels for the controls, standard and PRO, are both microprocessor operated but have different security levels. All the models of this series have open wire shelves, mounted on anti-tilt slides but offer the opportunity to fit drawers without any modification. The internal temperature is always guaranteed by the forced-air refrigeration system, with low noise cooling units. These are reliable and use environmentally friendly CFC-free gas.

Factory fitted options include lockable wheels and a chart recorder.

Accessories include pull-out drawers and extra shelves.



-45C Chest Freezers Also Available

Budget -40 Freezers

