



[Model GP75W](#)

FLUID WARMING CABINETS These are designed for use in theatres (urology, ENT, plastic surgery, cardiac, eye, day surgery) IVF, CAPD, dermatology, A&E, CT/Radiology and other areas where fluids, lotions or other media need to be kept warm or at body temperature. It is considered good practice to warm fluids prior to use as this helps maintain patient temperature, especially in critical areas such as Theatres, ITU and A&E.

DESCRIPTION

These cabinets will give long & trouble-free life in normal use. They are fitted with a digital or dial temperature controller & over-temperature safety cutout for product protection. Adjustable wire shelves allow free air movement around the stored bottles to guarantee better temperature

stability. Flat shelves are advisable when warming small bottles and bags. Special baskets are recommended when storing standard bags of fluid since they can be stored and accessed easier and additionally there is no risk of the outer wrapper slowing the movement of air round the cabinet or impairing the temperature stability or control.

The outer cabinets are made from sheet steel finished in easy-clean grey stove powder paint; stainless steel exterior cabinets are a standard option. Stainless steel chambers are more hygienic and durable, and offer excellent protection from rust and corrosion. Larger GP models are fitted with a handle (right) rather than a finger pull.

All models have an internal fan giving you optimum control and faster warm-up times. All models can be fitted with a glazed door panel and we can manufacture fluid warming cabinets to non-standard sizes and dimensions if required.



[Model GP30WB \(basic control\)](#)

Temperature Performance, all models:
 From 10c above ambient to +60c
 Stability typically +/- 0.2c
 Recovery after door opening is model dependant but is typically 1-3 minutes.



[Model GP150W fitted with glazed viewing panel](#)

FLUID WARMING CABINET RANGE [Combination fluid & blanket warming cabinets also available.](#)

Model Ref.	1 litre Bottles	3 litre Bags	No. Of Shelf / Basket	External (mm) Ht x W x D	Internal (mm) Ht x W x D	Weight Kg*
GP30W	Small bottles only		1 / -	570 x 415 x 422	220 x 315 x 315	20
GP40W	9	2	1 / -	670 x 415 x 422	320 x 315 x 315	25
GP50W	15	4	1 / -	615 x 590 x 422	240 x 490 x 315	30
GP52W	10	0	2 / -	845 x 590 x 300	555 x 490 x 193	35
GP75W	25	4	1 / 1	623 x 590 x 567	248 x 490 x 460	35
GP120W	30	4	1 / 1	735 x 590 x 659	360 x 490 x 552	45
GP130W	50	8	2 / 2	845 x 590 x 567	555 x 490 x 460	45
GP150W	60	12	2 / 2	845 x 590 x 659	555 x 490 x 552	50
GP210W	70	18	2 / 2	957 x 590 x 749	582 x 490 x 642	65
GP240W	90	24	3 / 3	1178 x 590 x 659	803 x 490 x 552	73
GP280W	98	32	2 / 2	941 x 750 x 769	566 x 650 x 662	78
GP350W	147	48	3 / 3	1104 x 750 x 769	814 x 650 x 662	95
GP355W	150	40	5 / 5	1720 x 590 x 659	1345 x 490 x 552	95

NOT suitable for warming blood products

Add 50mm to external width to allow for door opening.

* Empty