

# GENERAL PURPOSE LABORATORY INCUBATORS



## FEATURES:

Temperature ranges: from approx 5C above ambient; approx 10C above ambient for fan models.  
 Aluminium coated mild steel or stainless steel interior  
 Direct reading hydraulic controls with fluctuation of +/- 1.0C  
 Safety thermostat with indicator. Maximum temperature 100C

## OPTIONAL EXTRAS:

Fan assisted circulation (30 to 380 litre models)    Temperature calibration & mapping certificates  
 Digital control and readout (as picture) with fluctuation of +/- 0.1C    RS232 communications board  
 Programmable controllers and special timers (30 to 380 litre models)  
 Dial thermometer (30 to 350 litre models)  
 Inner glass door or glass outer door panel.    Door lock.    Left-hand hinged door  
 Operation on 110 volt supply

## DESIGN:

Exterior is sheet steel finished in easy clean stove powder paint  
 Incubator is fitted with shelf runners and removable chrome plated steel wire grid shelves  
 The door and chamber are fitted with silicone seals for excellent insulation against heat  
 The incubators are heated by Incoloy sheathed elements, 240v single phase supply

Model	Litres	Shelf Positions	No of Shelves	External mm			Internal mm			Max Load Height mm	Power watts
				H	x W	x D	H	x W	x D		
GP10	10	1	1	415	x 375	x 350	150	x 275	x 243	105	250
GP20	20	2	2	540	x 375	x 382	275	x 275	x 275	190	350
GP30	30	3	2	570	x 415	x 422	305	x 315	x 315	220	500
GP40	40	3	2	670	x 415	x 422	404	x 315	x 315	320	500
GP50	50	3	2	615	x 590	x 422	325	x 490	x 315	240	500
GP75	75	3	2	623	x 590	x 567	333	x 490	x 460	248	750
GP100	100	4	3	735	x 590	x 567	445	x 490	x 460	360	750
GP120	120	4	3	735	x 590	x 659	445	x 490	x 552	360	1000
GP150	150	5	3	845	x 590	x 659	555	x 490	x 552	470	1000
GP180	180	6	4	957	x 590	x 659	667	x 490	x 552	582	1250
GP210	210	6	4	957	x 590	x 749	667	x 490	x 642	582	1250
GP240	240	8	4	1178	x 590	x 659	888	x 490	x 552	803	1500
GP280	280	6	4	941	x 750	x 769	651	x 650	x 662	566	1500
GP350	350	7	4	1104	x 750	x 769	814	x 650	x 662	729	1500
GP380	380	8	4	1160	x 750	x 769	870	x 650	x 662	785	1500

Capacities up to 2500 litres and custom build – please call for alternative incubator leaflet.

**Standard Incubators**

Ringsstones Industrial Estate, Whaley Bridge, High Peak, SK23 7PD, UK  
 Tel: 01663 735494    Fax: 01663 732744    Web: www.gedscientific.co.uk



Drying Cabinets

Heavy Duty Ovens & Incubators

Cooled Incubators



Ovens

Special Build

