

GENERAL PURPOSE LABORATORY OVENS



FEATURES:

Temperature ranges: 40 to 300C (10 to 350 litre models). High temperature models 500C (50 - 250 litres)
 Aluminium coated mild steel or stainless steel interior
 Direct reading hydraulic controls with fluctuation of +/- 0.75C
 Safety thermostat with indicator

OPTIONAL EXTRAS:

Fan assisted circulation (30 to 350 litre models)
 Digital control and readout (as picture) with fluctuation of +/- 0.5C
 Programmable segment controllers
 Dial thermometer (30 to 350 litre models)
 Special timers (30 to 350 litre models)
 Operation on 110 volt supply

Incubator sterilizers
 Moisture extraction models

DESIGN:

Exterior is sheet steel finished in easy clean stove powder paint
 Oven is fitted with shelf runners and removable chrome plated steel wire grid shelves
 The door and chamber are fitted with silicone seals for excellent insulation against heat
 The ovens are heated by Incoloy sheathed elements, 240v single phase supply

Model	Litres	Shelf Positions	No of Shelves	External mm H x W x D	Internal mm H x W x D	Max Load Height mm	Power watts
GP10	10	1	1	415 x 375 x 350	150 x 275 x 243	105	350
GP20	20	2	2	540 x 375 x 382	275 x 275 x 275	190	500
GP30	30	3	2	570 x 415 x 422	305 x 315 x 315	220	750
GP40	40	3	2	670 x 415 x 422	404 x 315 x 315	320	1000
GP50	50	3	2	615 x 590 x 422	325 x 490 x 315	240	1000
GP75	75	3	2	623 x 590 x 567	333 x 490 x 460	248	1250
GP100	100	4	3	735 x 590 x 567	445 x 490 x 460	360	1250
GP120	120	4	3	735 x 590 x 659	445 x 490 x 552	360	1500
GP150	150	5	3	845 x 590 x 659	555 x 490 x 552	470	1500
GP180	180	6	4	957 x 590 x 659	667 x 490 x 552	582	1750
GP210	210	6	4	957 x 590 x 749	667 x 490 x 642	582	2000
GP240	240	8	4	1178 x 590 x 659	888 x 490 x 552	803	2250
GP280	280	6	4	941 x 750 x 769	651 x 650 x 662	566	2250
GP350	350	7	4	1104 x 750 x 769	814 x 650 x 662	729	2500

Capacities up to 2500 litres and custom build – please call for alternative incubator leaflet.



Drying Cabinets

Heavy Duty Ovens
& Incubators

Cooled Incubators



Incubators

Special
Build



Standard Ovens

Ringstones Industrial Estate, Whaley Bridge, High Peak, SK23 7PD, UK
 Tel: 01663 735494 Fax: 01663 732744 Web: www.qedscientific.co.uk

