

Trolley With Electric Height Adjustment

We supply stainless steel trolleys with electric height adjustment. These are ideal for loading & unloading washer disinfectors and autoclaves. They also have other uses within laboratories.

The trolley top can be fitted with any type of loading rails to make them compatible with all makes of washing and sterilizing equipment.

These trolleys are light compared with most others and this increases manoeuvrability and reduces manual handling issues. The trolley is fitted with four 80mm wheels with braking on the rear pair.

There are no hydraulics in the lift mechanism so there is no possibility of oil leaks; this can make them ideal for use in clean rooms. To make them acceptable for pharmaceutical use the 316 grade stainless steel is electro-polished.

In typical use the trolley can complete more than 50 operations before the battery needs re-charging (load dependant).

Supplied ready built.

Safe working load	190KG
Standard platform	545 x 870mm
Maximum height	900mm
Minimum height	330mm
Trolley weight	50KG
Re-charge time	12 hours

Other specifications - please call

QED Scientific Ltd.

Ringstones Industrial Estate
Whaley Bridge
High Peak SK23 7PD

Tel: 01663 735494
Fax: 01663 732744

Web: www.qedscientific.co.uk
E-mail: info@qedscientific.co.uk



Height Adjustable Tables & Desks



Electric or Manual Adjustment
Laboratory Tables & Desks
Cytology Desks
Microscope Desks
Clean Room Tables
Clean Room Trolleys
Autoclave Trolleys

Call for details of any product



C-Frame

Manual Height Adjustment

Correct working posture is critical where an operator rarely changes working position. Desks that the user can easily adjust are therefore to be recommended in order to comply with working height directives and reduce general fatigue.

For many applications manual adjustment of the table working height is acceptable.



I-Frame

Desks are available with either a cantilever or I-frame. We offer a wide selection of worktop dependant on the application. For general use a wood veneer finish is ideal, and we offer you a choice of wood grain finishes. For some uses wood is not acceptable and we offer alternatives including polypropylene, Trespa and stainless steel.



C-Frame with electric adjustment

Standard widths 1200, 1500, 1600, 1800 mm
Standard depths 750, 800 mm
Other sizes are available - please call.

Please refer to price list for full details.

Supplied flat packed for easy user assembly.

Wall Mounted (shown with sink)



In wet laboratory areas it is often recommended that as much equipment as possible is kept off the floor. For such areas we supply adjustable height tables that are fixed to a suitable load bearing wall so that the floor is kept completely clear. These tables can be fitted with a sink or other top as shown in the picture.



Radial corner desks



Scoop desks for microscopy

Electric Height Adjustment

Tables are manufactured from 304 or 316 grade stainless steel work top and foot assembly. This is important for many clean room applications.

The 4-leg model shown is designed specially for use with heavy items such as safety cabinets. The 2 leg version (below) suits most general purpose laboratory applications. Memory controls are fitted on the front as standard, but can be put elsewhere.

We are able to supply these to almost any size required and options include locking wheels, splash guard, fan/filter modules & power sockets.

Actuators are painted mild steel; stainless steel covers are an option at extra cost. Supplied ready-built for immediate use.

These desks are available in several widths and depths as well as a choice of work surface. Styles include rectangular, scoop and radial for corner positioning.

All models can be adjusted anywhere between 650 and 1200mm and lift up to 150KG. The front edges of the worktops are smoothed for comfort.

Worktops are available in a choice of veneers and include cherry, beech and maple. Oak is also an option in a limited number of sizes. These desks can also be fitted with alternative tops including Trespa, stainless steel and polypropylene.

For cytology screening there is a standard known as MDA02104, and these desks comply.

Supplied flat packed for easy user assembly.

Standard widths include 1200, 1600, 1800 & 2m. Standard depths include 800 & 1000mm. Other sizes are available so please call us. Wheels are a standard option on all models.

Please refer to price list for full details.