

These hot air sterilisers are dual purpose as they can be used as incubators as well as sterilizers, all having a high and low heat setting switch.



The exterior is constructed from sheet steel finished in easy clean stove powder paint. The interior can be a choice of aluminium coated mild steel (CLAD) or stainless steel (SS) and the shelves are chrome plated steel wire grids supported on fixed shelf runners.

The door and chamber are fitted with silicone seals for excellent insulation from heat and are easy to clean.

Incoloy sheathed heating elements are controlled by either a dial controller or digitally.

There are 14 models in the range with capacities from 50 to 350 litres.

Standard Specifications:

Temperature range 40C to 250C
CLAD interior
Direct reading hydraulic temperature control
Safety thermostat with warning indicator

Options Include:

Stainless steel chamber (SS)
Fan assisted air circulation (30 TO 35 L only)
Digital control and readout (DIG)
Dial thermometer readout
Programmable timers
Adjustable shelf ladders

Model	Litres	Shelf Positions	No of Shelves	External mm H x W x D	Internal mm H x W x D	Max Load Height mm	Power watts
GPIS10	10	1	1	415 x 375 x 350	150 x 275 x 243	105	250
GPIS20	20	2	2	540 x 375 x 382	275 x 275 x 275	190	350
GPIS30	30	3	2	570 x 415 x 422	305 x 315 x 315	220	700
GPIS40	40	3	2	670 x 415 x 422	404 x 315 x 315	320	700
GPIS50	50	3	2	615 x 590 x 422	325 x 490 x 315	240	1000
GPIS75	75	3	2	623 x 590 x 567	333 x 490 x 460	248	1000
GPIS100	100	4	3	735 x 590 x 567	445 x 490 x 460	360	1000
GPIS120	120	4	3	735 x 590 x 659	445 x 490 x 552	360	1500
GPIS150	150	5	3	845 x 590 x 659	555 x 490 x 552	470	1500
GPIS180	180	6	4	957 x 590 x 659	667 x 490 x 552	582	1500
GPIS210	210	6	4	957 x 590 x 749	667 x 490 x 642	582	2000
GPIS240	240	8	4	1178 x 590 x 659	888 x 490 x 552	803	2000
GPIS280	280	6	4	941 x 750 x 769	651 x 650 x 662	566	2000
GPIS350	350	7	4	1104 x 750 x 769	814 x 650 x 662	729	2500

Standard voltage 230v, 50/60Hz Other voltages to order