

**These are heavy duty commercial cabinets with a choice of sliding or hinged doors. These are equally at home in the pharmacy or laboratory.**



Model SD1140

**Standard Specifications**

Digital temperature display. \*  
 Internal circulation fan, with bottom refrigeration unit.  
 Adjustable shelves.  
 Fully automatic operation.  
 Electronic control, range +1°C to +12°C.  
 Easy clean white metal interior. \*  
 White/black exterior.  
 Auto-defrost (off-cycle). No drain required.  
 Hinged or sliding doors, double-glazed, self-closing. \*  
 CFC-free refrigerant (R134A).  
 Adjustable levelling feet.  
 Rated for use in high ambient temperatures to +43°C. \*  
 220/240v, 50/60Hz supply. 110v option – call.

Accessories include: extra shelves, door locks, and chart recorder.

**Alarm systems** – options include high & low alarm with delay. Audible & visual warnings with alarm mute & auto-reset. Please call.

**Sliding Doors**

Model SD890	675 litres	2 sliding doors, 8 shelves	2020h x 890w x 635d (mm)
Model SD1140	890 litres	2 sliding doors, 8 shelves	2020h x 1140w x 635d (mm)
Model SD1360	1075 litres	2 sliding doors, 8 shelves	2020h x 1360w x 635d (mm)
Model SD104	1040 litres	2 sliding doors, 10 shelves	1997h x 1110w x 820d (mm)
Model SD214	2040 litres	3 sliding doors, 15 shelves	1997h x 1970w x 820d (mm)

**Hinged Doors**

Model HD48CD	355 litres	1 hinged door, 4 shelves	1920h x 595w x 650d (mm)
Model HD690	525 litres	1 hinged door, 4 shelves	2020h x 690w x 635d (mm)
Model HD890	700 litres	2 hinged doors, 8 shelves	2020h x 890w x 635d (mm)
Model HD1140	920 litres	2 hinged doors, 8 shelves	2020h x 1140w x 635d (mm)
Model HD1360	1110 litres	2 hinged doors, 8 shelves	2020h x 1360w x 635d (mm)

\* HD48CD – plastic formed interior / hinged double glazed door / 30°C

**Chart Recorder**

The optional chart recorder is a 7-day disk type. It uses a standard long-life battery with a typical life of at least 1 year. The pens are replaceable fibre tipped. The recorder is complete with a pack of 100 charts, 125mm chart diameter, pens and see-through lockable door (2 keys are provided).  
 A mains powered model with battery backup is available on request.



*We also supply a wide range of solid door models for pharmacy and laboratory use. Please call.*

**Lab Glass Door Fridges**

Sycamore House, Eccles Road, Whaley Bridge, High Peak, SK23 7EW, UK  
 Tel: 01663 735494 Fax: 01663 732744 Web: www.qedscientific.co.uk

