




HDH style



Capacities between 200 and 1250 litres. Custom build.
Max temperature: 100C
Stainless steel chamber - option.
Digital microprocessor control.
Fan circulation as standard.
Heavy duty door hinges.
Safety thermostat with indicator.

OPTIONAL EXTRAS: Programmable controllers, Special timers
 Special shelves. 110volt supply.

DESIGN: Exterior is sheet steel finished in easy clean stove powder paint.
 Interior is fitted with shelf runners and removable chrome plated steel wire shelves.
 The incubators are heated by Incoloy sheathed elements, 240v single phase.

HDV style

RANGE:

Model	Litres	Shelves No. / Poss	External mm H x W x D	Internal mm H x W x D	Max Load Height	Power Rating
HDV200	200	3 / 4	1075 x 635 x 657	790 x 530 x 552	650	1000
HDV350	350	4 / 5	1165 x 795 x 758	880 x 690 x 650	740	1500
HDV425	425	4 / 7	1341 x 795 x 758	1056 x 690 x 650	916	2000
HDV500	500	5 / 8	1520 x 795 x 758	1235 x 690 x 650	970	3000
HDV625	625	4 / 6	1217 x 932 x 1008	932 x 827 x 900	792	2000
HDV750	750	5 / 7	1495 x 932 x 1008	1110 x 827 x 900	970	3000
HDV850	850	5 / 9	1520 x 932 x 1008	1233 x 827 x 900	1095	3000
HDV1000	1000	5 / 10	1847 x 932 x 1008	1462 x 827 x 900	1322	3000
HDH225	225	3 / 4	885 x 1065 x 550	600 x 960 x 445	460	2000
HDH350	350	3 / 4	885 x 1065 x 797	600 x 960 x 692	460	2500
HDH425	425	4 / 5	1003 x 1065 x 800	718 x 960 x 692	578	2000
HDH500	500	4 / 5	1121 x 1065 x 800	836 x 960 x 692	696	3000
HDH625	625	4 / 6	1217 x 1065 x 878	932 x 960 x 770	792	2000
HDH750	750	5 / 7	1495 x 1065 x 878	1110 x 960 x 770	970	4750
HDH850	850	4 / 5	1005 x 1645 x 954	720 x 1540 x 846	580	3000
HDH1000	1000	5 / 10	1847 x 1065 x 878	1462 x 960 x 770	1322	5400
HDH1220	1220	5 / 9	1475 x 1185 x 1138	1190 x 1080 x 1030	1050	4000
HDH1250	1250	5 / 10	1847 x 1252 x 908	1462 x 1147 x 800	1322	6000

As part of our policy of continuous product improvement, the specifications are subject to change without notice