

## MEDICAL & WARD REFRIGERATORS for storage of pharmaceuticals

Our 2007 range of medical refrigerators have been carefully designed to ensure that the internal storage temperature of between +2°C to +8°C is maintained at all times. This is achieved by a forced-air fan cooling system that ensures an even temperature distribution throughout the interior and rapid pull down & recovery after door openings.

Incorporated into the digital temperature display unit is an audio / visual alarm, which activates when a high or low temperature deviation is sensed. Integral to the unit is a min/max memory for continuous monitoring and manual daily recording of temperature.

In addition, these new medical refrigerator models have automatic defrosting to provide uninterrupted use without the need to defrost. Every medical refrigerator is fitted with a secure door lock for product security to prevent 'unauthorised' access to your valuable healthcare products.

To give you peace of mind this new range has been independently tested (UKAS Certified) to comply with WHO performance criteria (E3/PROC/3/2) for vaccine storage.

A 2 year UK warranty parts & labour guarantee covers all models in case of problems.

These refrigerators are also ideally suited for use on wards for storage of anything from pharmaceuticals to food products.

**Features:** +2°C to +8°C, Forced-Air Cooling, Digital Temperature Display, Audio / Visual Temperature Alarm, MAX/MIN Temperature Memory, Auto-Defrosting, Door Lock and keys, Meets W.H.O. performance criteria (E3/PROC/3/2), CFC free, Reversible Doors (Not SM40), Adjustable Shelving. Some models are available with the option of a glass door.

<u>RANGE</u>		<u>H x W x D mm</u>
<b>SM40</b>	42 litres	725 x 380 x 370
<b>SM60-G</b>	60 litres	650 x 450 x 430
<b>SM150</b>	150 litres	860 x 550 x 600
<b>SM180</b>	180 litres	900 x 600 x 600
<b>SM260</b>	260 litres	1225 x 600 x 600
<b>SM360</b>	360 litres	1595 x 600 x 600

### **OPTIONS** (factory fitted)

Glass door -G for model 60,180, 260 & 360  
Wall bracket for SM40  
Chart recorder  
Mains failure alarm  
Wheels



As part of our policy of continuous product improvement, the specifications are subject to change without notice.

**Medical & Ward Fridges**

Ringstones Industrial Estate, Whaley Bridge, High Peak, SK23 7PD, UK  
Tel: 01663 735494 Fax: 01663 732744 Web: [www.qedscientific.co.uk](http://www.qedscientific.co.uk)

**QED**  
Scientific