

# Telespro Rescue Covering

Protecting the casualty from heat loss and bacteria has significant benefits for further treatment. Simplified rescue process will also reduce the time spent in patient transport as well as improve the patient's safety.

The new Telespro Rescue Covering integrates the most essential properties of products for rescuing casualties and outperforms normal products used only for protection and patient transfer. It gives patients more intensive and effective protection throughout the entire rescuing process as well as makes the work of paramedics easier and more efficient.

The Telespro Rescue Covering is made of GORE™ Rescue Blanket laminate, which is particularly effective to protect patients against severe weather conditions, cold, moisture and wind, and bacteria penetration.

Telespro Rescue Covering for emergency care was

## Most common applications:

- Outdoor accidents
- Outbreak of sudden illness
- Covering not properly dressed patient
- Protecting casualties awaiting transport
- Rescue helicopters
- Rescuing near-drowning casualty
- Patients with burns
- Mountain rescue, e.g. at skiing centers


## Many functional properties

- Multiple zippers  
*Easy to put on: Freely opening and closing the covering when giving medical treatment to patient*
- Adjustable width  
*Flexible size that suits majority of the patients. Usable with neck braces, backboards, splints, etc*
- Reflectors  
*Ensure visibility and safety*
- Handles for carrying the patient  
*Easy to carry in constrained spaces e.g. staircase*
- Portable and foldable  
*Easy to move, carry and store*
- Pocket for using external heating pack  
*Facilitating the active warming of patient's body*

## Best quality materials

### GORE™ Rescue Blanket laminate

- Moisture and wind proof but breathable  
*Effective protection against cold, wind and moisture*
- Bacteria barrier  
*Safe for both patients and paramedics to use*
- Re-washable at 70° C  
*Hygienic*
- Highly durable for long time in use
- Soft textile surface  
*Relieving patients from discomfort*
- Insulating Underlay  
*Good thermal insulation that prevents patient from heat loss*



- 1 *Emergency service persons just leave the patient inside the covering at the hospital.*
- 2 *Emergency service persons will immediately get a ready to use covering from the stock of hospital.*
- 3 *Hospital will send the used covering to the hospital cleaning departments (hospital laundry)*
- 4 *Laundry will send the cleaned covering back to the hospital*

## Telespro's Maintenance Service

Keeping the ambulance and the equipment clean is of paramount importance for the safety of both the patient and the paramedics. With Telespro's Maintenance Service, the need to spend time cleaning up is reduced.

The unique characteristic of the distribution of our rescue coverings and services is that the rescue covers are leased for continuous use. This confers more benefits to the user; the initial investment is small and coverings are insured by Telespro and replaced if damaged.

The operation of the cleaning services creates a circulation of our covers between rescue services, hospitals and laundries.

We aim to improve:

**Safety** – Protecting both patients and rescue parties

Telespro's maintenance service guarantee the hygiene of rescue coverings:

- More effectively executing cleaning services and continuously provide customers with highly hygienic rescue coverings
- Clean coverings are always available for repeated use
- Marked and numbered covering allow Telespro to keep tracking and monitoring its maintenance services with high standard and specialized hospitals laundries

**Efficiency** – Reducing the time and variable cost of operation

Improved Efficiency of Rescue Operation shortens the time of rescue practice:

- Rescue departments will always get a clean covering immediately from the stock of hospital every time the patients are transferred, so that the time of operation is significantly reduced

- Rescue covering in use will remain in hospital and be sent to the laundry for cleaning and maintenance

**Continuous Protection** – Increasing the standard and effectiveness of rescue service

Ongoing protection is given to patients after arrival at the hospitals:

- Simplifying the patient transfer process
- Heat coverage continues without breaks even in hospital

Telespro's maintenance service model has been organized in Finland since spring 2004. We have received many favorable feedbacks from our customers in different field of rescue services, including ambulance emergency services; Finnish rescuehelicopter and coast guard rescue services; hospitals; and some private organizations, such as ski centers.

**UK Supplier:**

**QED Scientific Ltd.,** Wharf Road,  
Whaley Bridge, Derbyshire, SK23 7AD  
Tel: 0844 846 6435